



The Synergistic approach

Age-related anatomical changes are most effectively managed through the combination of science and technology. To achieve this goal, Neauvia research designed SMART COMBINATION THERAPY that combines different restorative therapies and grants clinically proven superior results.^{1,2}





2. Kolczewski et al., Hyaluronic Acid and Radiofrequency in Patients with Urogenital Atrophy and Vaginal Laxity, Pharmaceuticals 2022, 15, 1571, DOI: 10.3390/ph15121571



BASED DEVICES AND **COSMECEUTICALS**



NEAUVIA CVICCS

The best results with maximum safety

Neauvia provides a range of energy-based devices for the aesthetic, gynecological and surgical markets.

Neauvia devices are characterized by advanced electromedical engineering.

They are based on specifically designed and integrated technology, allowing optimal results in combination with Neauvia's PEG HA fillers and cosmeceuticals and ensuring maximum patient safety.^{1,2}



Kolczewski et al., Hyaluronic Acid and Radiofrequency in Patients with
Urogenital Atrophy and Vaginal Laxity, Pharmaceuticals 2022, 15, 1571, DOI:
10.3390/ph15121571



Focus on Zaffiro

Smart Hydro and Smart Infrared Tightening

Beautiful from every angle

Exfoliation of the skin

Non-invasive face lift

rejuvenation

Face and body Phototypes I to VI

ZAFFIRO IS THE ONLY MEDICAL DEVICE THAT COMBINES HYDRO-EXFOLIATION AND INFRARED THERMO-LIFTING TO BRIGHTEN AND TIGHTEN

- ► Waterpeeling to deliver precise hydroexfoliation, deep cleansing to the skin. It can be customized with 5 formulas: HA Skin Health, Acne formula, Soothing Skin formula, Skin Brightening and Hair nourishment.
- ► Infrared radiation to create gradual heating of the deeper layer of the dermis stimulating fibroblasts, to obtain a thickened, tauter and elastic skin.
- ▶ Up to 40% increase of fibroblasts and 25% increase on firmness 21 days post treatment.1

Zaffiro: Ila class Medical Device regulated health products which bears, under this regulation, the CE mark. Manufacturer: B&K, Domaniewska 37, Warsaw, Poland. Please carefully read the instructions in the leaflets. The use of this product requires the intervention of a healthcare professional 1. Data on file.



Part of Nlift **Smart Combination Therapy:** synergistic mid face protocol

Focus on Sectum

Multi-purpose RF Solution

Sectum

The shape of beauty



SECTUM HARNESS THE POWER OF RADIOFREQUENCY

How it works

Radiofrequency energy consists of alternating wavelength current that is delivered by one or more specially designed electrodes. As the energy travels across the tissue, the polarity or electrical charge of the skin tissue is changed. The friction of the movement of the energy current creates heat in the skin.

RF energy is emitted from one electrode and travels the path of least resistance or impedance in the skin and/or body, to a second electrode.

Sectum: Ilb class Medical Device regulated health product which bears, under this regulation, the CE mark. Manufacturer: B&K Domaniewska 37, Warsaw, Poland. Please carefully read the instructions in the leaflet. The use of these products requires the intervention of a healthcare professional. 1. Features and conditions of use in Sectum Product manual.

Aesthetic mode

Tighten and lift the face

Sculpt and shape the body



Part of Nboost Smart

Combination Therapy:

epidermis restoring

synergistic protocol



FDA

► High power - 120W maximum

- ► Low frequency 480kHz
- ► Smart Control technology for RF delivery optimization
- ► Temperature control from the handpieces
- ► Bipolar face applicator
- ► Multipolar body applicator

Gynecology mode

Tone and tighten the vulva skin and the vaginal canal



Part of Nrose Smart Combination Therapy: the 1st synergistic protocol dedicated to female intimate rejuvenation

- ► Bipolar applicator specially designed for gynecological treatments
- ► Non-invasive procedure
- ► Short processing time
- ► Temperature sensor integrated in the handpiece
- ► Smart Control technology for RF delivery optimization

Focus on Plasma

Portable Micro electrosurgery

Plasma IQ

Enhance every expression

Evelid Acne and skin laxity postoperative scars

lesion

PLASMA IQ IS AN ELECTRIC MICRO-SCALPEL DESIGNED FOR NON-INVASIVE SURGERY. IT USES PLASMA SUBLIMATION TO CAUSE A CONTROLLED AND VERY PRECISE SKIN LESION TO CREATE A THINNING AND RETRACTION OF THE SKIN TISSUE

- ▶ Ergonomic, cordless and portable
- ▶ 4 power levels
- ► Focused energy arc for precise treatment
- ► High quality electrodes for controlled energy delivery
- ▶ Low operating costs

Plasma IQ: IIb class Medical Device regulated health products which bears, under this regulation, the CE mark. Manufacturer: B&K, Domaniewska 37, Warsaw, Poland. Please carefully read the instructions in the leaflets. The use of this product requires the



Diode laser for hair reduction



Free to be free

Hair removal for all phototypes

EPILME IS EQUIPPED WITH
A LASER DIODE OF 808 NM
WAVELENGTH WHICH HAS
DEEP PENETRATION,
WHILE MAINTAINING A
HIGH SELECTIVE ABSORPTION
RATE BY MELANIN

- ▶ 4000W of power
- ► Handpiece guaranteed 100 million shots
- ► Constant cooling system at 1 °C
- ► Integrated pre-programming
- ▶ 3 modes of treatment

EpilMe: Ilb class Medical Device regulated health products which bears, under this regulation, the CE mark. Manufacturer: B&K, Domaniewska 37, Warsaw, Poland. Please carefully read the instructions in the leaflets. The use of this product requires the intervention of a healthcare professional.



Non-ablative skin resurfacing treatment



Unleash the best of your skin

Wrinkles

Skin Rejuvenation

Pigment treatments

Face and body



THE LATEST
NON-ABLATIVE
FRACTIONAL LASER
TECHNOLOGY
FOR SKIN
RESURFACING

- ▶ 1470 nm optimal wavelength
- ► Applicator window: 12x12 mm maximum
- ► Spot Spacing: 1-4 mm
- ► Effective in all skin types, including dark phototypes and sensitive skin
- ► Portable
- ► Low maintenance and technical assistance

LaserMe: Ilb class Medical Device regulated health products which bears, under this regulation, the CE mark. Manufacturer: B&K, Domaniewska 37, Warsaw, Poland. Please carefully read the instructions in the leaflets. The use of this product requires the intervention of a healthcare professional.



NEAUVIA PRECISION

Careful control and adjustment of energy delivery

- ► To optimize treatments for both results and patient comfort.
- ▶ Facilitating the use for medical staff.

NEAUVIA VERSATILITY

Comprehensive design to treat safely and effectively a variety of indications and areas.

NEAUVIA PRACTICALITY

Ease of use and cost-effectiveness that help providers offer world-class treatments at attractive prices.

In-House
Engineering & Concierge
Service



NEAUVIA®

- protocols, clinical results, techniques - on

www.neauviachannel.com





