




DIGITAL TOOTH DESIGN MADE BY CANDULOR

CANDULOR.

We have been designing
teeth since 1936.

Digital dentures for
removable restorations.





The 3Shape® Dental System is an innovative dental CAD design software which enables the digital design of dentures and the optimization of workflows. With the Candulor Denture SmileLibrary+ it now offers access to one of the largest tooth libraries, including set-ups for full dentures. The special feature is the predefined upper and lower teeth in occlusal relation. This way, the software accesses an idealized set-up in the database, enabling you to quickly find the optimal solution.

- 87 tooth molds: a comprehensive selection of anterior and posterior teeth in all available sizes
- 453 predefined set-up options of the angle class I
- Set-ups in the common occlusal schemes of Gerber, physiological, Gysi, lingualized and TiF
- Installation match in the multilayer structure of the TCR Multi disc

Set-up	Tooth molds	Candulor Denture SmileLibrary +	
		Upper Jaw anteriors	PhysioStar
			15
			PhysioSelect
			18
			PhysioSet
			22
		Lower Jaw anteriors	PhysioStar
			5
			PhysioSelect
			4
			PhysioSet
			8
		Upper and Lower Jaw Posteriors	Bonartic II
			3
			Bonartic
			5
			BonSelect
			4
			Condyloform II
			3
		Number of molds	Upper Jaw anteriors
			55
			Lower Jaw anteriors
			17
			Posteriors
			15
			Total
			87
		Occlusional concepts	Physiological
			✓
			Gerber
			✓
			Gysi
			✓
			TiF
			✓
			Lingualized
			✓
			Angle Class
			I
			Predefined Set-ups (Classe I)
			453

TOOTH MOLDS*

Upper Jaw Anteriors

PhysioStar	550, 552, 554, 556, 660, 662, 664, 666, 770, 772, 774, 776, 880, 882, 884
PhysioSelect	560, 562, 564, 566, 670, 672, 674, 676, 678, 680, 682, 780, 782, 784, 786, 788, 790, 792
PhysioSet	203, 204, 213, 224, 232, 233, 234, 262, 263, 264, 302, 303, 304, 332, 333, 342, 343, 353, 423, 424, 425, 452

Lower Jaw Anteriors


PhysioStar	990, 992, 994, 996, 998
PhysioSelect	60, 62, 64, 66
PhysioSet	52, 53, 54, 80, 81, 82, 83, 84

Upper & Lower Posteriors

Bonartic	02, 04, 05, 06, 08
Bonartic II	02, 04, 06
BonSelect	02, 04, 06, 08
Condyliform II	34, 36, 38

* Tooth moulds included in the Candulor Denture SmileLibrary+ offer

PROCESSING OPTIONS





Denture Bases

The denture base can be printed or milled, prefabricated teeth or milled tooth arches can be used.

Monoblocks

The base and teeth are made in one piece. It can be printed or milled as a try-in.



Tooth segments

Blocked tooth segments can be printed or milled. It is possible to design one order with several bridge segments.

PREDEFINED SET-UPS & OCCLUSION CONCEPTS

Anteriors	Posteriors	Gerber	Physio-logical	Gysi	Lingua-lized	TiF
PhysioSelect	BonSelect		4	4	4	4
	Bonartic		5	5		
PhysioStar	Bonartic II		3	3		
	Condylform II	3				

Posteriors	Gerber	Physiological	Gysi	Lingualized	TiF
BonSelect		✓	✓	✓	✓
Bonartic II		✓	✓		
Bonartic		✓	✓		
Condylform II	✓	✓		✓	



Wish to try?
We would be pleased
to advise you.

Candulor Dental GmbH


Am Riedergraben 6
D-78239 Rielasingen-Worblingen


T + 49 (0) 7731 79783-0

F + 49 (0) 7731 28917

E export@candulor.de

W candulor.com





768667/en/2025-01-01

CANDULOR. HIGH END ONLY.